



Euky Bear Steam Vaporiser

Special Features

- All night relief for colds and flu, bronchitis and croup
- Guaranteed quietest vaporiser on the market
- Revolutionary design which is safer and easier to use
- Works with any tap water
- Australian designed and owned
- Full 5 year warranty
- Suitable for children under 2 years



Quietest vaporiser on the market

The Euky Bear Steam Vaporiser is quiet, safer and easy to use. It is recommended for colds and flu, croup and bronchitis. It is guaranteed as the quietest vaporiser on the market. Using a safer water dispersible inhalant, it can be used Australia wide in all types of water. All night relief to breathe easier and sleep better, it is Australian designed with a full 5 year warranty.

Medical practitioners recommend the use of steam vaporisers. The steam moisturises mucous membranes, reduces irritation and clears mucus from the nose, bronchial tubes and lungs.

Steam Normalises Dry Air

In Winter, the relative humidity can drop to as low as 10%. To feel comfortable and stay healthy, we need to maintain a relative humidity of 40 – 50%. People with respiratory problems are better off when the humidity is increased to 60 - 80%.

Dry air causes the normally liquid thin protective coating of mucus on the linings of the airways to dry and thicken, becoming ineffective. When this happens, the delicate membranes of the nose, throat and larynx become irritated and highly susceptible to upper respiratory tract infections. The moist air from the vaporiser adds moisture to the air to lubricate and restore the use of these vital membranes.

Indications

- Helps relieve the symptoms of colds and flu, croup and bronchitis
- Helps relieve the severity of symptoms of colds and flu
- Helps relieve the symptoms of croup and bronchitis
- Helps provide humidity to assist the symptoms of cough including hoarseness
- Uses moist hot air to help thin phlegm and decrease viscosity and stickiness, helping phlegm to be expelled through coughing, thus reducing coughing
- Also helps create a healthier environment by adding moisture to the air

Special Features – Expanded

Quiet

The quietest steam vaporiser on the market. It is vital that patients, especially young children, sleep well and the whisper quiet operation enables patients to sleep soundly. Other vaporisers can be noisy and often disturb sleep.

Safer

Fully enclosed heating element. Some vaporisers have two electrodes in contact with the water. Whenever they are switched on the water is electrified and whenever water is electrified, there is the risk of accidental electrocution.

No medication cup. extralife Inhalant can be added directly to the water, greatly reducing the risk of accidental poisoning. Medication in the medication cup can attract young children. Children have dipped their fingers into the medication and then put their fingers in their mouths. This eliminates this problem.

No Special Water Required

Operates with any type of water suitable for drinking. It works perfectly in hard or soft water. The steam given off is absolutely pure. Any dissolved minerals in the water remain behind in the water container.



10-12 Hour Operation

Runs continuously for at least 12-14 hours irrespective of water hardness. Other vaporisers may promise long running but using hard water can run out of water in 4 hours or less, meaning it must be refilled in the night or the patient does not get the required treatment.

Starts Quickly

Starts steaming quickly, within three or four minutes with soft or hard water. Other vaporisers can take up to 30 minutes to get going.

Safety Night Light

Shows when the Euky Bear Steam Vaporiser is switched on. It is also useful to show exactly where it is in a dark bedroom.

Easy to Use

1. Pour some Eu-Clear Inhalant into the water container.
2. Add water to the required level. The water level is clearly seen through the side of the water container.
3. Lock the heating module into the water container.
4. Place the vaporiser into position.
5. Plug in and switch on power.

Cool and Efficient

Triple wall insulation of the heating module and concentric ridges in the base of the water container combine to prevent heat loss through convection from the core of the heating module.

Water in the water container remains cool during operation. The Euky Bear Steam Vaporiser is therefore safer to use and it is also more efficient because heat loss is reduced.

Easy to Maintain

The Euky Bear Steam Vaporiser is easy to maintain, clean and store after use. After each use, empty the water container, rinse the heating module and container and dry before storing. Dismantling, cleaning and re-assembly are unnecessary.

5 Year Warranty

The Euky Bear Steam Vaporiser carries a full five year warranty and is backed by a quick and efficient Service Centre.

Australian Designed

The Euky Bear Steam Vaporiser is Australian designed for Australian conditions. It is approved by Australian Electrical Authorities and it is registered as a medical device with the Therapeutic Goods Administration (TGA) Department of Community Services and Health Device Register.

Handy hints for getting the most out of your Euky Bear Steam Vaporiser:

- Use only water soluble inhalants, such as Eu-Clear Inhalant in your Vaporiser. Eu-Clear is water-soluble and specially formulated for your Euky Bear Steam Vaporiser.
- Use only Euky Bear Cleaning Tablets. Use of Euky Bear Cleaning Tablets helps protect the element and keep it free of mineral deposits, ensuring a longer life for your Vaporiser.
- Use fresh water each time the Vaporiser is used – do not “top up”. There is no need to add any other substances (ie. salt) to change the water quality as the Euky Bear unit operates in all kinds of water.
- Empty the water container and rinse and dry the heating module after each use to keep it clean and germ free.
- Always disconnect the cord from the power point before moving the unit (ie. to clean, store, fill or rinse).

Aust L: 11761.

Warning: If symptoms including coughing persist, consult your doctor.

Presentation and Availability

The Euky Bear Steam Vaporiser is available from pharmacies.